



Quality control standard metals for drinking water, low range, 500 mL

Artikel-Nr.: 2833749
CHF Preis: Kontakt
Verfügbar

Ready to use - no dilutions required. Mixed standard includes: Copper: 1 mg/l Cu, Iron: 0.3 mg/L Fe, Manganese: 0.1 mg/L Mn. 500 mL Bottle.

Technische Daten

Beschreibung:	Metals Drinking Water (Low Range) Mixed Parameter QC Standard
Konzentration:	Kupfer: 1,0 mg/L Cu
Plattform:	Lösung