



Filter membrane, EZ-Pak, 0.45 micron, 47 mm, white, gridded, 600/pk (millipore)

Artikel-Nr.: 2916800

CHF Preis: Kontakt

Kein Versanddatum angezeigt

Millipore EZ Pack Membranes, 0.45 µm pore size, 47mm

Dispenser provides superior convenience and assures sterility. Unique one press on the lever and a sterile membrane is presented (horizontally and inverted) for efficient and convenient assembly. The EZ Pak membranes come as a band of membranes that easily loads into the dispenser. Packaging materials are wound into the dispenser ensuring a clear operating area and quick change of bands. The membranes are Type HA, which signifies compliance with Standard Methods criteria.

Technische Daten

Durchmesser:	47 mm mm
Material:	Mixed-Cellulose Ester
Porengröße:	0,45 µm
Volumen / Packungsgröße:	600