



Membrane filters millipore S-pak type HA, 47 mm, 0.45 µm pore size, pk/600

Artikel-Nr.: 2936500
CHF Preis: Kontakt
Kein Versanddatum angezeigt

For membrane filtration in microbiological analysis of aqueous samples. "HA" signifies certified compliance with Standard Methods criteria with pore size of 0.45 µm. White, cellulose and mixed-ester composition. S-Pak (filter only) is sterile. Diameter: 4

Technische Daten

Durchmesser:	47 mm mm
Material:	Mixed-Cellulose Ester
Porengröße:	0,45 µm
Volumen / Packungsgröße:	600