



Einweg-BSB-Flasche, 300 mL

Artikel-Nr.: 2943101

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Single 300 mL disposable BOD bottle without stopper. Plastic bottle is EPA approved for BOD measurement.

Box of 100 bottles available under article number 2943100. Plastic stoppers (25 pack) available under article number 2943900. Glass stoppers (quantity of 1) available under article number 190901.

Technische Daten

Body Material:	Plastic
Elektroden-Typ:	Standard
Kapazität:	300 mL
Material:	PET
Volumen / Packungsgröße:	1 each