



Be Right™



Test strips total hardness, 0-425 mg/L, 50 tests

Artikel-Nr.: 2745250
CHF Preis: Kontakt
Verfügbar

Water hardness is caused almost entirely by calcium and magnesium ions. Other di- and trivalent metals have a similar effect, but usually are not present in high enough concentration in potable waters to cause problems. Hardness increases soap consumption in laundries and causes scale in boilers. Contains 50 tests.

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Technische Daten

Anzahl Bestimmungen:	50
Einheiten:	50 tests
Fußnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kleinste mögliche Schritte:	Steps: 0, 25, 50, 120, 250, 425 ppm
Kleinste Schrittweite:	Intervalos de: 0, 1.5, 3, 7, 15, 25 gpg
Messbereich:	0 - 425 ppm
Messbereich (2):	0 - 25 gpg
Parameter:	Hardness, Total As CaCO ₃
Plattform:	Teststreifen
Test-Kit Modell:	Test Strips
Versandgewicht:	0.2