




Be Right™



Rack, test tube, Polyethylene, 30 mm (O.D.), 21 holes

Artikel-Nr.: 2497904

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Twenty to ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Technische Daten

Durchmesser:	30 mm
Küvetten Anordnung:	3 x 7
Küvetten-Einsätze:	21
Volumen / Packungsgröße:	1 each