



m-Endo broth, dehydrated, 500 g

Artikel-Nr.:

1462300

Gefahrstoffe

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.



Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Technische Daten

Media Type:	m-Endo-Agar, dehydriert,
Region:	Global
Volumen / Packungsgröße:	500 g