



Be Right™



Bottle, sampling, presterilized, round, 50/pk

Artikel-Nr.: 2495050

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Polystyrene, disposable, transparent, non-fluorescent, break-resistant, presterilized bottles with caps. Bottles have 100 mL fill-to line; capacity of 175 mL. Bottles used with medium containing MUG or other fluorescing reagent. Work well with P-A Broth Ampules and PathoScreen(TM) Medium.

Technische Daten

Dechlorinating Agent?: No

Volumen / Packungsgröße: 50/pk