



TEST KIT, UNIV DIGITAL TITRATOR

Artikel-Nr.:

2270900

CHF Preis:

Kontakt

Verfügbar

Gefahrstoffe



Accurate, convenient titrations without the bulk, fragility, or waste of a conventional buret.

Universal Digital Titrator Kit with Digital Titrator, manual, 100-ml graduated cylinder, 125- and 250-ml Erlenmeyer flasks, demineralizer bottle, clippers, and rugged carrying case.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Technische Daten

Enthaltenes Zubehör:	Ja
Genauigkeit:	±1 % (Die Unsicherheit der Messwerte beträgt 1 Ziffer. Die meisten Proben erfordern mehr als 100 Ziffern.)
Gewicht:	0,132 kg
Inhalt:	Includes Digital Titrator, manual, 10-mL graduated cylinder, 125 - and 250-mL Erlenmeyer flasks, demineralizer bottle, clippers and a carrying case.
Lieferumfang:	Digital Titrator with accessories in carrying case

Inhalt

Includes Digital Titrator, manual, 10-mL graduated cylinder, 125 - and 250-mL Erlenmeyer flasks, demineralizer bottle, clippers and a carrying case.