



Wartungsset, Hoher Messbereich Phosphat-Analyser 5500sc/9611sc, 2/4 Kanal

Artikel-Nr.: 6788310

CHF Preis: Kontakt

Kein Versanddatum angezeigt

Maintenance kit includes: air filter, tubing, bottle caps with colors, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, user instructions

Technische Daten

Produkt-Kit: Hoher Messbereich Phosphat, 2/4-Kanal