



Be Right™



CDC565 Conductivity Cell, 4 poles

Artikel-Nr.:

E61M003

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Vorgeschlagene Ersatzprodukte:

- E61M010

Meets requirements of USP 24-NF19.

The right conductivity cell for your application.

It is important to choose a conductivity cell with the right construction and geometry for your particular application and working conditions. The Radiometer Analytical line offers a variety of conductivity cells to meet a large number of applications.

The right conductivity cell for your application.

It is important to choose a conductivity cell with the right construction and geometry for your particular application and working conditions. The Radiometer Analytical line offers a variety of conductivity cells to meet a large number of applications.