



Hardness, iron, and pH test kit, model HA-62A

Artikel-Nr.: 183701
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Model HA-62A

A multi-parameter test kit for hardness, iron & pH, this kit provides an easy-to-use, cost effective way to determine the quality of your water.

Popular tests at a great price in Hach's convenient and protective test kit carrying case.

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anwendung:	Water Quality
Anzahl Bestimmungen:	100
Messbereich:	6,6 - 8,4 pH-Einheiten
Modell:	HA-62A
Parameter:	Härte, Eisen und pH (6.6-8.4)
Plattform:	Color Disc