



Test kit, manganese color disc, model MN-PAN

Artikel-Nr.:

2350800

CHF Preis:

Kontakt

Kein Versanddatum angezeigt

Gefahrstoffe



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. Utilizing the simple color disc method, this kit measures manganese (Range: 0.05-0.7 mg/L) and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	50
Kleinste Schrittweite:	0.05
Koffertyp:	D
Messbereich:	0 - 0,70 mg/L Mn
Methode:	Color disc/PAN
Modell:	MN-PAN
Parameter:	Manganese
Plattform:	Color Disc
Standard-Messbereich:	0 - 0.70 mg/L Mn
Test-Kit Modell:	MN-PAN
Versandgewicht:	3