



Be Right™



## Test-Kit Kohlendioxid CA-23 1,25-25/5-40/5-100

Artikel-Nr.: 143601

CHF Preis: Kontakt

Verfügbar

Carbon dioxide test kit with drop count titration covers ranges of 1.25-25, 2-40, and 5-100 mg/L as CO<sub>2</sub>, with increments of 1.25, 2, and 5 respectively. Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

**Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.**

---

### Technische Daten

Anzahl Bestimmungen:	200
Fußnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kleinste Schrittweite:	1.25 2 5
Koffertyp:	D
Messbereich:	1.25 to 25 mg/L CO <sub>2</sub>
Messbereich (2):	2 to 40 mg/L CO <sub>2</sub>
Messbereich (3):	5 to 100 mg/L CO <sub>2</sub>
Methode:	Drop count titration/ Sodium hydroxide Phenolphthalein
Parameter:	Carbon Dioxide - As CO <sub>2</sub>
Plattform:	Tropfenzähltest
Test-Kit Modell:	CA-23
Versandgewicht:	1