



Dissolved oxygen test kit, model OX-2P

Artikel-Nr.: **146900**
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Low dissolved oxygen levels in surface water systems usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. This test kit utilizes drop count titration and contains 100 tests.

Method / Chemistry: Drop count t

Utilizes drop count titration - a simple, economical titration method

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.2
Koffertyp:	D
Messbereich:	0.2 - 4.0 mg/L & 1 - 20 mg/L O ₂
Methode:	Drop count titration/ Modified Winkler
Parameter:	Dissolved Oxygen - As DO
Plattform:	Tropfenzähl-Titration
Test-Kit Modell:	OX-2P
Versandgewicht:	3