



Be Right™



Sulfite test kit, model SU-5

Artikel-Nr.:

148002

CHF Preis:

Kontakt

Verfügbar

Gefahrstoffe



Utilizing drop count titration/iodometric to determine sulfite, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen: 100

Kleinste Schrittweite: 1

10

Koffertyp:

D

Messbereich:

1.0 - 20 mg/L & 10 - 200 mg/L SO₃²⁻

Methode: Drop count,

titration Iodometric

Parameter: Sulfite - As SO₃²⁻

Plattform: Tropfenzähl-Titration

Test-Kit Modell: SU-5

Versandgewicht: 2