



Test kit, ozone, model OZ-2

Artikel-Nr.: 2064400
CHF Preis: Kontakt
Lieferung innerhalb von 2 Wochen

Gefahrstoffe



Hach's most economical test kit for ozone measurement utilizes the simple, low-cost color disc method.

Hach's most economical test kit for ozone measurement

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.05
Koffertyp:	B
Messbereich:	0 - 2.2 mg/L O ₃
Methode:	Color disc/DPD
Modell:	OZ-2
Parameter:	Ozone
Plattform:	Color Disc
Standard-Messbereich:	0 - 2.2 mg/L O ₃
Test-Kit Modell:	OZ-2