



Be Right™



Flask, Filtering, Polypropylene, 1000 mL

Artikel-Nr.: 2089953

CHF Preis: Kontakt

Kein Versanddatum angezeigt

Polypropylene filtering flask with side-arm at 45° angle prevents tipping. Barbed side-arm accepts 1/4-3/8" I.D. tubing for application of vacuum. See page 380 for tubing; see page 375 for stopper.

Technische Daten

Einheiten:	1 Stück
Gehäusematerial:	Polypropylen
Kapazität:	1000 mL