



Bottle, Storage, Wide Mouth, Write-On, 500 mL, 12/pk

Artikel-Nr.: 2930003

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Large dedicated labeling areas help to standardize the way bottle contents are identified. A separate area is provided for name and a larger area allows for detail of contents or safety information and includes a defined space for date information. Bottles are molded of rugged HDPE in 500 and 1000 mL sizes and are available as narrow or wide mouth with leak proof polypropylene cap.

Technische Daten

Body Material:	High Density Plastic
Kapazität:	500 mL
Material:	HDPE
Package size UOM:	12pk
Volumen / Packungsgröße:	12 /bag