



CEL 890 Professionelles portables Trinkwasserlabor

Artikel-Nr.:

2688500

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Bestehend aus:

- Tragbarem Colorimeter zur chem. Wasseruntersuchung (DR890),
- Tragbarem pH/Leitfähigkeit-Messgerät mit pH-Elektrode und 4-Pol LF-Elektrode (HQ40D.99.12201),
- Digitaltitrator
- Reagenzien-Erstausrüstung
- Zubehör
- Illustrierte Bedienungsanleitungen

Komplett im Tragekoffer.

Contains durable instruments, premeasured reagents impervious to environmental contamination, and apparatus and labware

Offers on-site results in an easy-to-use, rugged, portable lab

Rugged kit is designed to withstand harsh field conditions

Hard-sided durable cases

Instruments have option of running on batteries, so tests can be run immediately without concern for sample preservation and handling

versatile for lab or field use

Real-time decisions

In the field, rapid testing and decision making is essential. With this kit, there is no waiting for hours or days for a lab result

Reliable results

Whether used in your mobile lab or at the field site, this CEL890 kit yields consistent results time after time

Technische Daten

Anwendung:

Potable Water

Inhalt:

The Advanced Portable Laboratory includes DR/890 Colorimeter, HQ40d multi-meter, PHC201 pH electrode (1 meter cable), CDC401 conductivity probe (1 meter cable), Digital Titrator, illustrated instrument and procedures manuals, reagent sets, and all required apparatus in a rugged carrying case.

Inhalt

The Advanced Portable Laboratory includes DR/890 Colorimeter, HQ40d multi-meter, PHC201 pH electrode (1 meter cable), CDC401conductivity probe (1 meter cable), Digital Titrator, illustrated instrument and procedures manuals, reagent sets, and all required apparatus in a rugged carrying case.