



## ORP Platinum combination electrode

Artikel-Nr.:

5193700

### VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

A good performer in dirtier, more difficult samples where redox measurements are commonly made. This platinum-type electrode with temperature sensor features a reference junction that resists clogging and is easy to clean. Includes two gel cartridges and

### Accurate

You don't just need an answer, but the right one, and fast. Optimized stabilization algorithms in Hach meter platforms eliminate operator guesswork and reduce inaccurate measurements. Many electrodes work specifically to deliver accurate results and t

### Easy to Use

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

### Rugged and Durable

Rely on Hach's durable measuring systems designed to withstand the rigors of outdoor testing and frequent use. Choose from the industry's toughest probes or non-glass electrodes for rugged, waterproof systems that deliver day after day.

### Flexible

Need to measure in both the lab and field? Hach offers efficiency with reliable, high-performance systems that easily move with you from the benchtop directly to the testing site. Measure the parameters you need today with the assurance you can expand them tomorrow.

---

## Technische Daten

Diaphragma:	Open
Durchmesser:	12.3 mm
Elektroden-Typ:	Laboratory
Elektrolytlösung:	2546902
Gewährleistung:	6 Monate
Kabellänge:	0.91 m
Länge:	241 mm
Material:	Sensor Body: PMP
Messbereich:	± 2000 mV
Sensortyp:	Platinum annular ring
Spezielles Ausstattungsmerkmal:	Platinum series
Temperaturbereich:	N/A N/A N/A

---

## Erforderliches Zubehör

- Kaliumchlorid Referenz-Elektrolyt Gel-Kartusche (Item 2546902)