



pH combination electrode, refillable

Artikel-Nr.: 5194011

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Standard combination pH probe with a refillable reference and built-in temperature sensor. Contains fixed cable with BNC and 3.5 mm phone connectors compatible Hach One series meters.

Accurate

You don't just need an answer, but the right one, and fast.

Easy to Use

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

Rugged and Durable

Rely on Hach's durable measuring systems designed to withstand the rigors of outdoor testing and frequent use. Choose from the industry's toughest probes or non-glass electrodes for rugged, waterproof systems that deliver day after day.

Flexible

Need to measure in both the lab and field? Hach offers efficiency with reliable, high-performance systems that easily move with you from the benchtop directly to the testing site. Measure the parameters you need today with the assurance you can expand them tomorrow.

Technische Daten

Connector (text):	BNC and 3.5 mm Phone Connectors
Diaphragma:	Glass frit
Durchmesser:	12.1 mm
Elektroden-Typ:	Laboratory
Elektrolytlösung:	5194100
Gewährleistung:	6 Monate
Kabellänge:	1 m
Länge:	152.5 mm
Material:	Sensor Body: Epoxy
Messbereich:	0 - 14 pH
Sensortyp:	Glass
Temperaturbereich:	N/A N/A N/A
Thermistor:	# (ATC)