



Sodium chloride, 10,000 $\mu\text{S}/\text{cm}$, 500 mL

Artikel-Nr.: 2972249

CHF Preis: Kontakt

Kein Versanddatum angezeigt

Conductivity standard solutions for laboratory and field analyses

Für Leitfähigkeit-Elektroden, ebenfalls geeignet für die Einstellung von Leitfähigkeit-Apparaturen, 10000 $\mu\text{S}/\text{cm}$ ($\mu\text{mhos}/\text{cm}$)

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Technische Daten

Beschreibung:	Conductivity calibration solution
Gebindeeinstellungen:	Bottle
Konzentration:	10000 $\mu\text{S}/\text{cm}$ at 25°C $\pm$ 100 $\mu\text{S}/\text{cm}$ at 25°C
Lagerbedingungen:	10 - 25 °C (protect from light)
Parameter:	Leitfähigkeit
Plattform:	Lösung
Volumen / Packungsgröße:	500 mL