



Filter paper, glass fiber, extra thick, pore size 1 µm, diameter 47 mm, 100/pk

Artikel-Nr.: 2842100

CHF Preis: Kontakt

Kein Versanddatum angezeigt

For filtration of heavily contaminated samples. High particulate-holding capacity makes discs efficient depth filters and allows for filtration of large volumes of solutions. Pall Life Sciences.

Technische Daten

Durchmesser:	47 mm mm
Material:	Glass Microfiber
Porengröße:	1 µm
Volumen / Packungsgröße:	100