



Beaker, Polypropylene, Low Form, 250 mL, 6/pk

Artikel-Nr.: 108076

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Technische Daten

Durchmesser:	83.6 mm with spout
Einheiten:	6 Stück
Gehäusematerial:	Polypropylene
Kapazität:	250 mL
Material:	Polypropylen, niedrige Form