



Test kit, nitrogen, nitrate color disc, model NI-14

Artikel-Nr.: **1416100**
CHF Preis: Kontakt
Lieferung innerhalb von 2 Wochen

Gefahrstoffe



Utilizing the simple color disc method, this kit measures Nitrogen, Nitrate as NO₃ and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	50
Kleinste Schrittweite:	0.02
	0.2
Koffertyp:	B
Messbereich:	Dual Range
Messbereich (2):	0 - 10 mg/L NO ₃ ⁻ -N
Methode:	Color disc/ Cadmium reduction
Modell:	NI-14
Parameter:	Nitrogen, Nitrate
Plattform:	Color Disc
Standard-Messbereich:	0 - 1 mg/L NO ₃ ⁻ -N; 0 - 10 mg/L NO ₃ ⁻ -N
Test-Kit Modell:	NI-14**
Versandgewicht:	1