



## Test kit, copper (free) test color disc, model CU-5

Artikel-Nr.:

1421300

CHF Preis:

Kontakt

Kein Versanddatum angezeigt

Gefahrstoffe



Copper may exist in natural waters and effluents as a soluble salt or as suspended solids. A small amount is essential for plants and animals, and concentrations exceeding 0.1 mg/L are useful for controlling algae and plankton growth. Quantities ranging from 0.02 to 0.1 mg/L are toxic for some fish, so its use for treating fish disease requires careful monitoring. This kit is USEPA-accepted for waste water and potable water and approved by some states for potable water under specific conditions (Check with your regulatory agency). Range: 0.1-5 mg/L; contains 100 tests.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Technische Daten

Anzahl Bestimmungen:	100 F
Kleinste Schrittweite:	0.1
Koffertyp:	B
Messbereich:	0 - 4 mg/L Cu
Methode:	Color-disc/Bicinchoninate
Modell:	CU-5
Parameter:	Copper (free)
Plattform:	Color Disc
Standard-Messbereich:	0 - 4 mg/L Cu
Test-Kit Modell:	CU-5
Versandgewicht:	1