



Iron and manganese color disc test kit, model IR-20

Artikel-Nr.: 146300
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



When customers complain of stained and dingy laundry, this kit lets you do quick, on-the-spot tests for iron and manganese to see where the problem lies.

Method / Chemistry: Color disc colorimetry. Test: Iron/ Manganese. Range (mg/L): 0-5/ 0-3. Smallest I

Simple, fast, color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anwendung:	Water Quality
Kleinste mögliche Schritte:	0.2
Messbereich:	0 - 3,0
Messbereich (2):	0 - 4,0 mg/L Fe
Modell:	IR-20
Parameter:	Eisen und Mangan
Plattform:	Color Disc