



Hydrazine color disc test kit, model HY-2

Artikel-Nr.: 184900
 CHF Preis: Kontakt
 Kein Versanddatum angezeigt

Gefahrstoffe



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method with a range of 0.02-1 mg/L, this kit contains 300 tests.

Measures Hydrazine, a common oxygen scavenger in boiler water, using the simple color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	300
Kleinste Schrittweite:	0
Koffertyp:	D
Messbereich:	0 - 1 mg/L N ₂ H ₄
Methode:	Color disc/ p-Dimethylamino- benzaldehyde
Modell:	HY-2
Parameter:	Detergents
Plattform:	Color Disc
Standard-Messbereich:	0 - 1 mg/L N ₂ H ₄
Test-Kit Modell:	HY-2
Versandgewicht:	2,5