



Bottle, Sampling, Low Density, Polyethylene, 30 mL, 12/pk

Artikel-Nr.: 2087050
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Low-density polyethylene (LDPE) bottles feature wide-mouth opening for easy filling and leak-proof polypropylene closure. Great choice for grab samples, shipping liquids, soil sampling, and storing dry chemicals.

Technische Daten

Body Material:	Low Density Plastic
Kapazität:	30 mL
Material:	LDPE
Package size UOM:	12pk
Volumen / Packungsgröße:	12 /pk