



## Centrifuge Tube, Conical Bottom, 10/pk

Artikel-Nr.: 2115839

CHF Preis: Kontakt

Kein Versanddatum angezeigt

Centrifuge tube with conical bottom. Polycarbonate with molded-in 0.1mL graduations. For low-speed centrifugation up to 6000 x g. Autoclavable. The 15 mL size is also used with Settlometer Kit (P/N 21157-00). Supplied in packages of 10, without closures. Autoclavable. Nalge.

---

### Technische Daten

|                          |        |
|--------------------------|--------|
| Kapazität:               | 15 mL  |
| Länge:                   | 120 mm |
| Volumen / Packungsgröße: | 10 /pk |