



## Bag, sterile, WhirlPak, 207 mL, 25/pk

Artikel-Nr.: 2233125  
CHF Preis: Kontakt  
Kein Versanddatum angezeigt

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection. "Write-on" surface accepts waterproof ink. 24 oz bag perfect for incubating petri dishes in water bath applications.

---

### Technische Daten

Kapazität: 207 mL  
Volumen / Packungsgröße: 25