



## Test kit, water conditioning, model SO-1

Artikel-Nr.: 225300

### VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Hardness/Soap "Demo" for Sales Aid in Water Conditioning. Use drop-count methods to demonstrate the relationship of hardness to soap consumption. Note: While it provides a quick hardness/soap check, it is less accurate than the hardness titration supplied by Model 15-E (PN 148201) and Model AD-16E (PN 148410) water conditioning test kits.

**Uses drop-count methods to demonstrate the relationship of hardness to soap consumption**

**Provides an easy, economical, and fast hardness/soap check**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Technische Daten

Messbereich:	1 - 20 gpg
Parameter:	Soap Demonstration
Plattform:	Drop Count Titration