



Test kit, chlorine (total), model CN-65

Artikel-Nr.: 225401
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	100
Kleinste mögliche Schritte:	1
Kleinste Schrittweite:	0.2
	1
Koffertyp:	D
Messbereich:	0,2 - 4,0 mg/L & 1 - 20 mg/L
Methode:	Drop count titration/ Thiosulfate
Parameter:	Chlorine, mid and high range - As total Cl ₂
Plattform:	Tropfenzähl-Titration
Test-Kit Modell:	CN-65
Versandgewicht:	2