



Reagent set, chromate, digital titrator

Artikel-Nr.: 2272400

Gefahrstoffe

VERALTETER ARTIKEL
Dieser Artikel ist nicht mehr verfügbar.



For determination of chromate by sodium thiosulfate Digital Titrator titration. Hach Method 8211. Range: 20-400 mg/L CrO₄. Reagent Set includes: Dissolved Oxygen 3 Reagent Powder Pillows, Potassium Iodide Powder Pillows, Sodium Thiosulfate Cartridge (0.2068 N) and Starch Indicator Solution. Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Technische Daten

Anzahl Bestimmungen:	40 - 100
Messbereich:	20 - 400 mg/L CrO ₄ ²⁻
Methode:	8211