



Test kit, hydrogenperoxide, model HYP-1

Artikel-Nr.: 2291700
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Model HYP-1 Hydrogen Peroxide Test Kit (as H_2O_2) utilizes drop count titration and contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.2
	1
Koffertyp:	D
Messbereich:	0,2 - 2,0 mg/L & 1 - 10 mg/L
Methode:	Drop count titration/ Thiosulfate
Parameter:	Hydrogen Peroxide - As H_2O_2
Plattform:	Tropfenzähl-Titration
Test-Kit Modell:	HYP-1
Versandgewicht:	2,5