



Reagent solution pillows, silver 2, pk/50

Artikel-Nr.: 2293666

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Gefahrstoffe



Reagent for determination of Silver by the Cadion 2B Colorimetric Method. Hach Method 8120. Range: 0.02 - 0.70 mg/L Ag. Sample Size: 50 mL. Pack of 50 Pillows. Also available in reagent set 2296600.

Technische Daten

Anzahl Bestimmungen: 50

Methode: Colorimetric