



Reagent, FerroZine® iron, 25 g

Artikel-Nr.: 230424
CHF Preis: Kontakt
Kein Versanddatum angezeigt

[3-(2-Pyridyl)-5,6-bis(4-phenylsulfonic acid)-1,2,4-triazine, monosodium salt, monohydrate] CAS No. 69898-45-9. One of the most sensitive iron reagents known $\text{N}_4\text{C}_{20}\text{H}_{13}\text{S}_2\text{O}_6\text{Na} \cdot \text{H}_2\text{O}$, FW = 510.47. Assay: >95%

Technische Daten

Methode: FerroZine
Plattform: Behälter