



Be Right™



Butterfield's buffered phosphate

Artikel-Nr.: 2319110

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Butterfields Buffered Phosphate, 450 mL bottle. pk/10.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Technische Daten

Einheiten:	10 St.
Volumen:	450 mL
Volumen / Packungsgröße:	10/pk