



Filter, acrocap, pore size 0.2 µm

Artikel-Nr.: 2323810

CHF Preis: Kontakt

Kein Versanddatum angezeigt

For sterilization or clarification of up to three liters of culture media or other aqueous solutions. Ideal for sterilization of solutions that tend to foam when filtered under vacuum. Inlet accepts 6.4 mm (0.25") tubing or most male luer fittings.

Technische Daten

Durchmesser:	n.z.
Filter Typ:	Filter, acrocap
Packing:	n.z.
Porengröße:	0,2 µm