



Be Right™



Test kit, boiler feed & scale, model BSC-1, drop count titration

Artikel-Nr.:

2350300

CHF Preis:

Kontakt

Kein Versanddatum angezeigt

Gefahrstoffe



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. An economical test kit for monitoring scale formation in small industrial boilers, the BSC-1 Boiler Feed and Scale Test Kit uses drop count titration methods to determine levels of alkalinity, chloride, and hardness.

Technische Daten

Anwendung:	Boiler Water
Modell:	BSC-1
Plattform:	Tropfenzähltest