



Test kit, molybdate, model MO-LR

Artikel-Nr.: 2359300
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method

Effectively monitors corrosion inhibitors

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.25
Koffertyp:	D
Messbereich:	Single Range
Methode:	Color disc/ Ternary complex
Modell:	MO-LR
Parameter:	Molybdate
Plattform:	Color Disc
Standard-Messbereich:	0 - 3.0 mg/L Mo
Test-Kit Modell:	MO-LR
Versandgewicht:	2