



Buffer solution pillows, phosphate, sterile 50

Artikel-Nr.: 2375866

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Phosphate buffer solution. AODC Procedure for Total Bacteria. pk/50 pillows.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Technische Daten

Volumen / Packungsgröße: 50/pk