



Filter, acrodisc, pore size 0.2 µm

Artikel-Nr.: 2415600
CHF Preis: Kontakt
Kein Versanddatum angezeigt

For filtration of dilute biological fluids. Sterilized by gamma irradiation to eliminate potential contamination by EtO residuals. Inlet/Outlet connections: 13/25 mm: Female Luer-Lok inlet, standard male slip luer outlet. Pore size: 0.2 µm. Diameter: 25mm.

Technische Daten

Durchmesser:	NA
Filter Typ:	Filter, acrodisc
Material:	Polysulfone
Packing:	NA
Porengröße:	0.2 µm