



Probenbehälter, 45 mL

Artikel-Nr.: 2552114
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Plastic, leak-proof. Corning Snap Seal™ offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with patented sealing mechanism so sample won't leak.

Technische Daten

Durchmesser:	30 mm (approx)
Kapazität:	45 mL
Volumen / Packungsgröße:	600 /pk