



Thermometer, Enviro-Safe, Non-Mercury, -10 - 260 °C

Artikel-Nr.: 2635701

CHF Preis: Kontakt

Kein Versanddatum angezeigt

Enviro-Safe non-Mercury thermometers feature an exclusive biodegradable liquid and dye. Easy to read, aesthetically distinctive due to dark green liquid and permanently fused markings. Each thermometer, individually calibrated against NIST traceable equipment, identified by unique serial number for ISO

9000 or other quality systems.

Technische Daten

Länge:	405 mm
Temperaturbereich:	-10 °C - 260 °C
Thermometer-Einteilung:	1 °C
Thermometer-Flüssigkeit:	Non-Mercury
Volumen / Packungsgröße:	1 each