



Cylinder, graduated, mixing, certified, 25 mL

Artikel-Nr.: 2636340
CHF Preis: Kontakt
Kein Versanddatum angezeigt

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Technische Daten

Gehäusematerial:	Glass
Kapazität:	25 mL
Messmethode:	To Contain
Typ:	Single Metric, Mixing
Volumen / Packungsgröße:	1 each
Zertifizierung:	Klasse A