



Vacuum filter microfunnel, pk/50.

Artikel-Nr.: 2831500

CHF Preis: Kontakt

Kein Versanddatum angezeigt

Feature a completely assembled filtration unit: Funnel, filter, nutrient pad and petri dish. Filter: GN-6 Metrical with 0.45 µm pore size and 47 mm diameter; white, gridded surface. Funnel: Polypropylene Petri Dish: polystyrene lid and polypropylene base Nutrient Pad: Cellulose. Fits 2831500, 2846100 and 2846150 manifolds. Sterile and individually bagged.

Technische Daten

Durchmesser:	47 mm
Material:	Mixed-Cellulose Ester
Porengröße:	0,45 µm
Volumen / Packungsgröße:	50