



Quality control standard inorganics for drinking water, 500 mL

Artikel-Nr.: 2833049

CHF Preis: Kontakt
Lieferung innerhalb von 2 Wochen

Inorganics Quality Control Standard for Drinking Water. Ready to use. No dilution required. Mixed standard includes: Fluoride: 1 mg/L F^- , Nitrate: 2 mg/L NO_3^- -N, Phosphate: 2 mg/L PO_4^{3-} , Sulfate: 50 mg/L SO_4^{2-} . 500 mL Bottle.

Technische Daten

Beschreibung:	Drinking Water Inorganics, Mixed Parameter QC Standard
Konzentration:	Fluorid 1 mg/L F^-
Plattform:	Lösung