



Quality control standard hardness, low range, 500 mL

Artikel-Nr.:	2833449
CHF Preis:	Kontakt
Verfügbar	

Ready to use - no dilutions required. Mixed standard includes: Total Hardness: 100 mg/L as CaCO₃. Calcium Hardness: 50 mg/L as CaCO₃.

Technische Daten

Beschreibung:	Low Range Hardness, Mixed Parameter QC Standard
Konzentration:	Calciumhärte 50 mg/L CaCO ₃
Plattform:	Lösung